# Eye Pose Tool – Full Baseline Technical Design Document (Bilingual)

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- Ophthalmologists, orthoptists, clinical researchers.
- 眼科医生、斜弱视与视功能医生、临床研究人员。

### 1.3 Core Value / 核心价值

- Automatic cropping with a fixed 2:1 aspect ratio (眉毛 鼻梁，上下；脸廓 耳廓，左右)。
- Automatic gaze direction classification into 9 directions (UL, U, UR, L, C, R, DL, D, DR) regardless of input order.
- Automatic placement into a fixed 3 × 3 grid (B-area), with optional manual drag correction.
- High-quality JPEG export (quality=100, subsampling=0) suitable for medical records.

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---

### 2.1 High-level Workflow / 总体流程

1. User selects 1 – 9 photos in any order    A-area shows raw thumbnails.
2. Flutter calls Python to process all images:
  - Face & eye detection
  - Eye pose normalization (roll correction + centering)
  - 2:1 cropping
  - Gaze direction estimation (9 directions)

3. Flutter auto-fills B-area  $3 \times 3$  grid based on directions.
4. Doctor may drag images in B-area to correct misplacements.
5. User taps “ Save Photos ” to export a final  $3 \times 3$  mosaic to the gallery.

## 2.2 A-area (Raw Input) / A 区 ( 原始输入 )

- $3 \times 3$  grid showing up to 9 original photos.
- The order has \*\*no meaning\*\* for directions.
- A-area is purely for input preview.

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## 2.4 Gaze Direction Estimation / 视线方向判定

- For each eye, compute dx, dy:  
 $dx = (pupil\_x - eye\_center\_x) / eye\_width$   
 $dy = (pupil\_y - eye\_center\_y) / eye\_height$
- Use thresholds  $T_x=0.25$ ,  $T_y=0.25$  to classify Left/Center/Right and Up/Center/Down.
- Combine to 9 directions: upLeft, up, upRight, left, center, right, downLeft, down, downRight.

## 2.5 B-area Auto Layout / B 区自动布局

- Fixed layout:  
Row1: upLeft | up | upRight  
Row2: left | center | right  
Row3: downLeft | down | downRight
- Each GazeResult is mapped to one cell according to its direction.
- If multiple images share same direction: keep the one with higher confidence, others go to a candidate list (optional).

## 2.6 Manual Drag Correction / 拖拽纠错

- Long-press to drag a cell in B-area; dropping onto another cell swaps the two.
- Manual operations override automatic classification for final export.

## 2.7 Export / 导出

- Combine B-area current layout into one  $3 \times 3$  high-quality JPEG.
- Save to Android gallery; the exported mosaic must visually match B-area exactly.

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-----

## 3.1 Overall Layout / 整体布局

- Single scrollable page (SingleChildScrollView).
- Layout from top to bottom:

A-area (3 × 3 raw grid)

B-area (3 × 3 processed grid, draggable)

C-area (3 buttons: Select / Process / Save)

### 3.2 A-area / A 区

- 3 × 3 non-scrollable grid, 9 cells.
- Shows original photos or placeholders with “ + ” icon.
- Does not indicate gaze directions.

### 3.3 B-area / B 区

- 3 × 3 grid, fixed semantics per position:  
(0,0)=upLeft, (0,1)=up, (0,2)=upRight,  
(1,0)=left, (1,1)=center, (1,2)=right,  
(2,0)=downLeft, (2,1)=down, (2,2)=downRight.
- Displays processed 2:1 cropped images.
- Supports long-press drag-and-drop to swap cells.

### 3.4 C-area Buttons / C 区按钮

- Select Photos: open system picker, allow 1 – 9 selection, refill A-area and clear B-area.
- Process Photos: send current A-area images to Python, fill B-area with results; disabled if no images.
- Save Photos: export B-area mosaic; disabled if B-area has no processed images.

### 3.5 State Management / 状态管理

- Recommended: Riverpod.
- Core states:
  - rawPhotos: list of original paths.
  - processedGrid: List length 9.
  - processingState: idle / processing / done.
  - errorMessage: last error to show.

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### 4.1 Responsibilities / 职责

- Read input images.
- Detect face & landmarks (Mediapipe FaceMesh/Iris).

- Normalize head pose (roll correction; eye-eye line horizontal).
- Center nose line & eyes for cropping.
- Compute gaze direction dx/dy.
- Crop to 2:1 region (眉毛 鼻梁, 左右脸廓).
- Save processed images.
- Return JSON with direction, confidence, paths.

#### 4.2 Face & Iris Detection / 人脸与虹膜检测

- Use MediaPipe FaceMesh with refine\_landmarks=True.
- Key regions (indices may vary by version):
- Left eye contour
- Right eye contour
- Iris points for each eye
- Nose bridge top & nose tip

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#### 4.4 Gaze Vector Calculation / 视线偏移计算

- For each eye, compute pupil center and eye center.
- $\text{dx} = (\text{pupil\_x} - \text{eye\_center\_x}) / \text{eye\_width}$
- $\text{dy} = (\text{pupil\_y} - \text{eye\_center\_y}) / \text{eye\_height}$

#### 4.5 Direction Classification / 方向分类

- Thresholds: Tx=0.25, Ty=0.25.

- 水平 :  $dx < -Tx$  Left,  $dx > Tx$  Right,  $|dx| \leq Tx$  Center。
- 垂直 :  $dy < -Ty$  Up,  $dy > Ty$  Down,  $|dy| \leq Ty$  Center。
- Combine into 9 directions: upLeft, up, upRight, left, center, right, downLeft, down, downRight.

#### 4.6 Confidence Score / 置信度

- Based on:
  - Iris detection quality.
  - Consistency between left and right eye dx/dy.
  - Magnitude of deviation beyond thresholds.
- Range: 0.0 ~ 1.0.

#### 4.7 JSON Response / JSON 返回结构

Example:

```
{
  "results": [
    {
      "original_path": "/input/p1.jpg",
      "cropped_path": "/output/p1_crop.jpg",
      "direction": "upLeft",
      "confidence": 0.92,
      "rotation_angle": -3.2
    }
  ],
  "errors": []
}
```

#### 4.8 Error Codes / 错误码

- "NoFaceDetected"
- "IrisNotFound"
- "InvalidGazeVector"
- "LowConfidence"
- "CropFailed"
- "RotateFailed"
- "PythonCrash"
- "FileMissing"
- "UnsupportedFormat"

### 5.1 Communication / 通信方式

- Use `serious_python` as the bridge between Flutter and embedded Python.
- Flutter calls:  
`await SeriousPython.invoke("process_images", argsJson)`

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---

#### 7.1 Flutter Modules

- ui/: home\_page.dart, grid\_a\_area.dart, grid\_b\_area.dart, button\_panel.dart, image\_viewer.dart
- models/: gaze\_result.dart, eye\_direction.dart, process\_result.dart
- providers/: raw\_photo\_provider.dart, grid\_provider.dart, python\_state\_provider.dart
- services/: python\_service.dart, export\_service.dart

#### 7.2 A-area Module

- Displays raw photos.
- No semantics about directions.

#### 7.3 B-area Module

- Displays processed & auto-sorted photos.
- Supports drag-and-drop.

#### 7.4 Drag Module

- Implemented via LongPressDraggable + DragTarget.
- Swaps GazeResult entries in grid list.

#### 7.5 Python Module

- processor.py as main logic.
- face\_detector.py, iris\_detector.py, gaze\_vector.py, direction\_classifier.py, cropper.py.

#### 7.6 Export Service

- Reads B-area grid.
- Asks Python (or Dart) to compose final 3 × 3 JPEG.
- Saves to gallery.

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- Show “ 请先选择照片 / Please select photos first. ”
- Disable Process button.

### 8.2 NoFaceDetected

- Skip that image, keep B-area cell empty.
- Show message including filename.

### 8.3 IrisNotFound, InvalidGazeVector, LowConfidence

- Do not auto-place into B-area.
- Optionally keep as candidate; doctor can correct by drag.

### 8.4 Multiple Images Same Direction

- Keep the higher confidence in main grid.
- Lower confidence ones go to candidates.

### 8.5 CropFailed, RotateFailed

- Show human-readable error; skip that image.

### 8.6 PythonCrash

- Show “ 处理失败，请重试 / Processing failed, please retry. ”
- Optionally restart Python engine.

### 8.7 Empty Cells in Export

- Still export full  $3 \times 3$  mosaic; empty cells use gray background.

### 8.8 Very Large Input

- Automatically downscale before processing.

### 8.9 OnlyOneEyeDetected, EyesClosed

- Mark as low confidence or error; require manual correction.

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-----

### 9.1 Performance Targets

- 9 images processed    3 seconds on a flagship device.
- Python engine init    1.2 seconds.
- Memory usage    400 MB peak.

### 9.2 Python Optimization

- Load Mediapipe model once.
- Batch process all images.
- Downscale very large images.

### 9.3 Flutter Optimization

- Use Riverpod to minimize rebuild.
- Use lightweight drag widgets.

### 9.4 Device Compatibility

- Android 10 – 15.
- ARM64-v8a as primary ABI.
- No dependency on Google Play Services.

### 9.5 Offline Requirement

- Entire system must run fully offline.
- No cloud calls, no external APIs.

---

## Chapter 10. Future Extensions / 第 10 章 未来扩展

---

### 10.1 AI-based Gaze Classifier

- Replace or augment geometric dx/dy with a small CNN model.

### 10.2 Quality Assessment

- Score photo quality: blur, darkness, occlusion.

### 10.3 Multi-patient Management

- Store multiple mosaics with patient IDs & timestamps.

### 10.4 Annotation Tools

- Allow drawing arrows, notes on each image.

### 10.5 Video Mode

- Capture a short video, auto-select best frames for each direction.

---

End of Baseline Document / 基线文档结束

---